



INFORMATION CARD


TCG | TEMPERATURE CONTROL GATE

www.hs-solutions.eu

01

Using the device

MEASUREMENT POSITION




TCG | Temperature control gate


www.hs-solutions.eu

RECOMMENDATIONS

01 You need to take off your glasses




02 You need to uncover your forehead




CORRECT POSITION +/- 0,25 °C

Parallel to the screen

**NOT RECOMMENDED POSITIONS +/- 0,5 °C**

Leaning away from the screen



Leaning towards the screen

